



Properties	Rapid Setting -1	Rapid Setting -2	Slow Setting -1	Slow Setting -2
Residue on 600 micron IS Sieve, % Max	0.05	0.05	0.05	0.05
Viscosity @ 50 ° C(Saybolt Viscometer), Sec	20-100	100-300	20-100	30-150
Storage Stability after 24 Hours, %	Less than 2	Less than 1	Less than 2	Less than 2
Binder residue by evaporation, % Min.	60	67	50	60
Coagulation at low temperature	Nil	Nil	Nil	Nil
Miscibility with water	No Coagulation	No Coagulation	Immiscible	No Coagulation